Germa Rescue Lite Full Body Splint

Rel.300420

Distributor:

Ferno Washington Inc.
70 Weil Way
Wilmington, OH 45177 - USA

Phone: +1 937 382 1451
Fax: +1 937 382 6569

E-mail: info@ferno.com
www.ferno.com

Manufacturer:

AB Germa
Industrigatan 54-56
SE-291 36 Kristianstad - Sweden

Phone: +46 (0)44 12 30 30
Fax: +46 (0)44 10 31 79

E-mail: info@germa.se
www.germa.se

*Read this manual
carefully and keep it
for future reference*

Germa Rescue Lite Full Body Splint

English - Directions for use_____	Page	4
Deutsch - Bedienungsanleitung_____	Seite	8
Français - Mode d'emploi_____	Page	12
Español - Instrucciones de uso_____	Página	16
Svenska - Bruksanvisning_____	Sida	20
Italiano - Istruzioni per l'uso_____	Pagina	24
Dansk - Brugsanvisning_____	Side	28
Nederlands - Gebruiksaanwijzing_____	Bladzijde	32
Suomi - Käyttöohje_____	Sivu	36
Norsk - Brukerveiledning_____	Side	40

ENGLISH - DIRECTIONS FOR USE

1. Introduction

The Germa Rescue Lite is a full body splint. The splint is designed to be placed on top of a Basket stretcher/trolley, as a patient mattress. The full body splint is principally designed for the immobilisation of patient where leg and/or spinal injuries are suspected. Deep tissue damage can also be avoided by using full body splint.

The full body splint is equipped with 3 handles on both sides. It is also equipped with four patient restraint straps, each with Fix-Lock, quick-release plastic buckles.

The full body splint should be operated with a minimum of 2 medical trained operators when carrying the full body splint for shorter distance. Additional people may be required for heavier patients and when carrying longer distance. Additional carrying belts can also be used to aid manual handling.

The full body splint can be dragged on the ground using the carry handles or by connecting a carrying belt or rope to the full body splint, making sure to avoid rough surfaces or sharp edges on the ground along the route of movement.

The Germa Rescue Lite is a full body splint consists of an outer shell made from airtight material, inner chambers system which contain granules. When removing the air from full body splint using the supplied vacuum pump, the shell becomes rigid and can be shaped around the patient as the air in the splint is removed.

When the air has been fully evacuated from the Germa Rescue Lite full body splint, it is rigid and it is not possible to mould it any more. The mattress wraps around the patient offering an optimal support and stabilization, without exerting any pressure on the patient.

The Germa Rescue Lite full body splint can also be used to support, float aid of a patient in water. (See precautions section).

Recommendation: Usable life of 5 years

WARNING

These directions for use may be updated without further notice. Copies of the current version are available from the manufacturer.

2. Precautions

- Never mould the Germa Rescue Lite full body splint on the top of the head or under the foot soles as it might cause undesirable pressure against the spine.
 - On a regular basis check that the vacuum pump is functioning correctly, and physically test it by deploying the full body splint by removing air from it until the shell and mattress becomes hard.
 - Avoid handling the splint next to sharp objects. If the splint is being dragged along the ground, be aware of obstacles on the ground, such as sharp stones, sticks, glass, sharp edges etc.
 - Check pump, full body splint and carrying straps for any visible damage or wear and tear. Repair any damage as required before re-using the full body splint. Change any defect parts immediately. Ensure that all parts are present and operating, including the valve cap, patient restraints and buckles and head strap
 - Avoid handling the splint next to open fire.
 - On a current basis check that the splint has no visible cracks.
 - The Rescue Lite full body splint is unstable in water. You must never leave a patient unattended in water at any time.
 - When transporting in an aircraft, please note that the atmospheric pressure in the cabin may drop when reaching higher altitudes and the vac splint will therefore become less rigid. In this situation, make sure to evacuate more air out from the full body splint.
- Own modifications on product not allowed – Your medical specialist to be responsible for any modification.

When decided to take the product out of service and destroy – make sure to follow your local instructions or contact us for additional information.

3. Description

The Germa Rescue Lite consists of an outer shell made of hard-wearing PVC coated fabric and a unique granulate-filled inner bag. When the air is evacuated from the vacuum splint the individual granulates pack up together into a tight mass, making the splint rigid. Before and while the air is evacuated the splint may be moulded and positioned according to the shape of the patient, to give best immobilisation. An optimum support of the patient is obtained when all air has been evacuated. This can be determined by the fact that it becomes more difficult to pump as less air remains in the vac splint, to the point where the pump cannot be moved any more once all the air has been removed.

The full body splint is equipped with 3 handles on both sides. It is also equipped with four colour coded patient restraint belts with Fix Lock buckles. When not in use, store the belts and fasten them by use of the buckles in the gap between the row of handles and the mattress. When the Rescue Lite full body splint is not being used, and is not required to be used as a splint on top of a Basket stretcher / trolley, the vacuum splint can be folded and packed to store in the supplied protective bag together with pump and repair kit, taking up less space.

4. Operating the Rescue Lite Vacuum Splint

1. Put the vacuum splint on an even surface possible near the patient, trying to avoid any sharp objects that may be on the ground. Check that the valve is closed by turning clockwise until fully closed, if needed then twist anti-clockwise $\frac{1}{2}$ a turn to open the valve. Note that the cap should not be removed from the valve at any time during operation, to avoid it getting lost.

2. Position and adjust the splint on the ground so that there is an even distribution of granules along all areas of the vac splint, and that it is flat ready to take the patient, with restraint straps unbuckled and not underneath the vac splint

3. Check that the valve is closed by turning it full clockwise until it feels tight, and then connects the pump.

4. Place the injured patient in the medically desired position on the full body splint, utilising standard adopted manual handling procedures, utilising supplementary equipment and slide sheets to aid in safe transfer of the patient.

5. Make sure to support stabilisation of the patient head and neck. Mould the splint around the patient body to give effective immobilisation, such that it becomes stable and comfortable to the patient. Mould up some of the granular material in the base of the mattress between the legs of the patient. Ensure that the head flaps are positioned effectively around the patient's head and ensure the head strap is attached

6. Connect the patient restraint straps according to the colour code, from legs to chest, and tighten them one by one. X Chest belts are red and light blue. Cross them over the chest.

WARNING

Make sure that the vac splint is not pressing on top of the patient's head or on the soles of their feet

7. Evacuate the air using the vacuum pump and continue moulding the splint around the patient's body as more air is removed and as the splint becomes more rigid. Tighten the restraint straps further as the vac splint becomes stiffer to give better moulding around the patient. Check that there is no additional pressure to the top of head or under the feet. Check that the unit is stable, give good immobilization and is comfortable for the patient.

8. The patient is now ready to be transported, use the handles for patient transport and handle with care.

4.1 Use of the Rescue Lite in combination with Basket stretcher

The Rescue Lite can be used as a soft Stretcher/Transfer Mattress, when the air is not removed from the unit.

1. Ensure that the splint valve is opened to allow air to flow into the splint.
2. Place the patient on the full body splint, fold the side sections of the splint over patient and use the patient restraints to secure the patient on the splint
3. Place the Rescue Lite with the patient in desired position of the Basket stretcher and secure the patient on the stretcher prior to start transfer process.

5. Cleaning and Disinfecting

After each use the vac splint should be wiped with a moist cloth and then dried. If dirty the splint should be washed with soapy water and rinsed with clean water. Allow the splint to dry completely before storage or usage.

If the mattress needs disinfecting first clean the splint carefully and let it dry. Wipe the complete mattress with a solution containing 70% alcohol or most common disinfectant can be used e.g. the disinfectant e.g. called LYSETOL AF.

6. Maintenance and Repair

At frequent regular intervals the splint, valve and pump must be checked for any visible damage.

At the same time the splint must be tested for function and the suction operation of the suction pump must be checked according to Section 7.

Leaks should be checked for by visual inspection or by slightly filling of the mattress with air, covering it with soapy water all over the surfaces and the welds, then watch for air bubbles that may appear, which will indicate a leak. When filling the mattress with air to check for leaks, ensure that the mattress is not filled so much such that it creates a pressure within the splint, which may then cause the mattress edge weld to split.

If any leak is found on the splint, first clean the area around the leak and allow it to dry.

Find a repair patch of a suitable size from the repair kit and then spread the supplied glue all over the patch and all around the area of the splint where the patch will be applied.

Wait for the glue to become "dust dry" (approx. 2-3 minutes) and then place the patch on the area of glue and leakage.

Apply pressure to the whole area of the patch; making sure good contact has been established in the whole area and that any air bubbles are removed.

The repair will have full strength after 24 hours, but can be used within a few minutes in case of emergency

Note: Patient restraint straps are detachable, washable and replaceable, washing in laundry bag, see label at the foot end restraint. If blood gets on the straps, soak in cold water first before washing.

7. Testing the Vacuum Splint and Pump

Straighten out the splint on the floor, close the valve by turning it fully clockwise and connect the pump.

Evacuate the air until the splint becomes rigid and leave it for about 2 hours. Check to ensure that the splint is still as rigid. Any loss of rigidity may indicate a leak in the vacuum splint or through the valve. Check both for leaks and correct functioning of the valve cap.

See Section 6 Maintenance and Repair or check the valve for leakage from the vacsplint.

If OK, try to insert a new valve and repeat the test again.

If the pump does not seem to be removing air from the vacuum splint, check the pump function by removing it from the vacuum splint, place a finger over the orifice at the end of the pump tubing and pump a few times. A negative

Germa Rescue Lite Full Body Splint

pressure should be created in the tubing and should be felt at your fingertip.

If the pump does not function do not try to repair it, return it for service to your local dealer or to a Germa service centre.

8. Storage

After each use, check to ensure that the full body splint is still operating correctly and made ready for its next use. Open the valve by turning it fully clockwise to fully close then 1/2-1 turn anti-clockwise to open and put the splint onto an even surface. Make sure the granules have been distributed evenly, close the valve by turning clockwise and then fold the mattress in its initial position.

The vacuum splint can be stored on top of the stretcher, or folded in the storage bag, such that it will fit in the designated storage area.

9. Specifications

The Rescue Lite Full Body Splint is in conformity with Council Directive 93/42/EEC concerning Medical Devices MR Safe.

Dimensions

Length x Width x Height	see product information below
Weight	see product information below
Heat resistance	70°
Cold resistance	-30°

Capacity for foot pump

The pump has capacity to empty the mattress according to EN 1865.

Materials

Part	Material
Top and bottom material	Polyester fibres PVC coated (flame resistant)
Extra reinforcement material	Polyester fibres PVC coated (flame resistant)
Innerbag	Polyester
Wood	Birch

Product

Cat. no

271 513 000

271 513 001

Product name

Germa Rescue Lite Full Body Splint

Germa Rescue Lite Full Body Splint Ready2go

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Report no. :

(To be assigned
in by Manufacturer)

Please print or type or fill in by MS Word, all sections a. – i. Fax or e-mail to your supplier.
Please use repair order when reporting repairs not associated with complaints.

For claims please specify :
Purchase Date

Description of product

Description of Incident/Problem: 1) What, 2) How?

When? In connection with: 1) transport, 2) unpacking, 3) preparing, 4) test, 5) cleaning, 6) normal use, 7) training,
8) disassembly, 9) assembly, 10) storage, 11) other:

Germa Rescue Lite Full Body Splint

Info on the product, prior to problem,

Storage: 1) stock, 2) site of use, 3) packed, 4) unpacked, 5) softpack, 6) other:

Used cleaning and disinfecting agents/Sterilizing method:

Pattern of use: 1) stand-by for emergency, 2) routine use, 3) stand by, 4) other:

Environment: 1) Indoor, 2) Outdoor, 3) Heat, 4) Cold, 5) Humidity, 6) Sunlight, 7) Gases, 8) other:

Which **actions** have been: a) taken, b) planned by the reporter or user at the site of incident?

1) examination of the product on site, 2) replacement of the product, 3) other:

Suggested actions by the distributor/complainant:

Returning of samples/faulty product(s),

1) 1) Already shipped/Will be shipped to _____,

2) 2) Pending request from manufacturer,

3) 3) Product/parts not available:

Date

We acknowledge the receipt of this report by manufacturer and will start processing it immediately. Please use the assigned complaint report number for further reference.

Your contact person in this matter is: _____

Name:

Direct Tel.

Direct Fax.

Init/Date:

Rel. 300420

Distributor:

Ferno Washington Inc.
70 Weil Way
Wilmington, OH 45177 - USA

Phone: +1 937 382 1451
Fax: +1 937 382 6569

E-mail: info@ferno.com
www.ferno.com

Manufacturer:

AB Germa
Industrigatan 54-56
SE-291 36 Kristianstad - Sweden

Phone: +46 (0)44 12 30 30
Fax: +46 (0)44 10 31 79

E-mail: info@germa.se
www.germa.se 12442400301, 2015-01-20, version 3