Stretcher Instructions

EN [Product Name] Stretcher [Model] AZ-CFS, AZFST-01, AZ-FST-02, AZ-FST-03, AZ-FST-04, AZ-SCP, AZ-SCA, AZ-BSK, AZ-BSK2, AZ-STR, AZ-LSTR-01, AZ-LSTR-02, AZ-LSTR-03

[Intended Use] For the transport of the transferred patients.

[Application] Camping adventure, emergency rescue, disaster prevention reserve, etc.

[Usage environment] Outdoor camping, adventure, emergency, mountain climbing, camping.

[Usage Method] Spread the stretcher flat on the ground, put the patient flat on the stretcher, and lift the handles at both ends of the stretcher.

[Cleaning and Disinfection] 1. Clean all surfaces of the stretcher with a cloth, warmwater and soap or a mild detergent soap.

2. If necessary, to remove stains, use a soft brush or a lightsolvent, paying attention to labels. 3. Dry with a cloth. Make sure that all the components of the device are completely dry before reusing them. Do not wash with high pressure water. It is advisable to wear suitable gloves (e.g. household disposable gloves) during disinfection. **[Storage]** Avoid high temperature & moisture **[Shelf life]** 5 years

 **Attention** Always check whether each part of the folding stretcher is flexible and firm, safe and reliable. If problems are found, it should be repaired and replaced in time.

NO [Produktnavn] Båre [Modell] AZ-CFS, AZFST-01, AZ-FST-02, AZ-FST-03, AZ-FST-04, AZ-SCP, AZ-SCA, AZ-BSK, AZ-BSK2, AZ-STR, AZ-LSTR-01, AZ-LSTR-02, AZ-LSTR-03

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[Labels, Packing Logo Design]

Symbol	Introductions	Symbol	Introductions
	CE Symbol		AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	medical device		AUTHORIZED REPRESENTATIVE IN THE SWITZERLAND
	Manufacturer Name Address		Manufacture Date
	Batch Code		Use-by date
	Warnings and Precautions		Consult instructions for use
	Catalog Number		Do not use if package is damaged
	Unique Device Identifier		Importer

Anping Guardian Medical Equipment Co., Ltd.
500 meters West of Nan Su Village, Anping County, Hengshui City, 053600, Hebei, P.R. China.
Tel: +86 13932089888 www.med-resq.com
info@tjresq.com

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MedNet EC-REP GmbH
MedNet SWISS GmbH