



# FUNERAL & MORTUARY

# INTRODUCTION

## *The Ferno Story – A Legacy of Innovation and Ergonomics*

**Since 1955, FERNO has been dedicated to delivering high-quality Funeral and Mortuary solutions, supporting every aspect of your work—from first call to final presentation.**

### *Shaping the Future of Patient Transport*

The evolution of modern ambulance cots is deeply rooted in the need for better ergonomics for caregivers and greater comfort and dignity for patients. In the early 1900s, stretchers were basic two-pole designs, often used to transport the deceased—commonly by funeral homes.

By the 1920s and '30s, wheeled stretchers began to reduce physical strain. The breakthrough came in 1952 with the first elevating cot, making bed-level transfers possible and marking the start of ergonomic innovation.

But the true turning point came in 1955, when Dick Ferneau founded Ferno Manufacturing Company. Just a few years later, in 1958, Ferno introduced the game-changing One-Man Mortuary Cot—a revolutionary design that allowed a single operator to safely and efficiently handle removals. This innovation sparked a new era of ergonomic, user-friendly stretchers and set the standard for the entire industry.

Today, FERNO is the global leader in the design and manufacture of equipment for the mortuary, emergency, and rescue sectors—driven by the same commitment to safety, dignity, and innovation.

Our full line of equipment, from lifting, operating, and dressing tables to durable, easy-to-clean church trucks, is built to make your job easier and your service more professional. Every FERNO product is designed to the highest industry standards—because we believe in helping you serve with confidence and care.

***If you have any questions or need assistance, we're here to help.***

***Let's work together to find the best solutions for your ergonomic needs.***

***Contact us today!***



## VIPER-F – Power Mortuary Stretcher with Loading System

The VIPER-F is a high-performance electric mortuary stretcher designed for modern, professional funeral services. With a load capacity of 320 kg and a lifting height up to 90 cm, it offers ergonomic, intuitive operation and exceptional maneuverability.

The VIPER-F includes a wide lift-off Stretcher Model 11BS with adapters to also carry a coffin. One product that handles both transports of bodies on a lift-off stretchers or in a coffin.

- Electrohydraulic lifting system – smooth raising and lowering at the push of a button.
- Ergonomic and user-friendly – intuitive +/- controls and soft-start/stop function.
- Durable and easy to maneuver – aluminum frame with 150 mm swivel wheels and full rubber soft tires.
- Automatic loading system – vehicle integration allows effortless loading/unloading, even on uneven surfaces.
- Crash-tested vehicle mount included – with telescopic guides and auto-lock for secure fixation.
- Integrated tablet/smartphone connectivity for setting height parameters and monitoring system status.
- Onboard battery system – automatically recharges when docked in the vehicle.

### Adapter V 161 S – Multi-Function Coffin & Stretcher Adapter

- Secure stretcher lock – prevents vertical lifting during transport.
- Dual mode operation:
- With brackets folded down: transports FERNO heavy-duty stretcher 11 BS.
- With brackets upright: securely holds a coffin.
- Quick-mount system – easily attaches to the frame and locks into place.
- End-stop plate (M 3369) – adds extra safety during stretcher transport.

### Vehicle Loading System (Standard Equipment)

- Easy and smooth insertion/removal with ball-bearing telescopic guides.
- Mechanical centering allows angled entry.
- Manual override system for loading/unloading without power.
- Battery charging starts automatically when docked.

### 11BS Lift-Off Stretcher

- 10 cm wider than regular transfer stretchers
- Extra lightweight yet loadable up to 250 kg
- Designed specifically for the transport of very heavy individuals



<b>Art. nr</b>	F1030
<b>Heigh</b>	38-133 cm
<b>Length</b>	160-200 cm
<b>Load Capacity</b>	320 kg (400 kg)
<b>Weight</b>	75 kg (66kg)

# MORTUARY STRETCHERS

## Ferno Mondial Funeral Transporter

The Ferno Mondial Funeral Transporter is designed for use with our mortuary stretchers and is perfectly suited for the dignified transport of coffins and deceased individuals. Thanks to its innovative SMART-LOAD™ technology, there's no need to manually operate a lever when loading or unloading the roll-in system into or out of a vehicle—the release mechanism is fully automatic. With the adapter included, our mortuary stretchers can be securely and easily attached. The adapter features two crossbars with stretcher supports and a fixed stop plate, ensuring a stable and reliable connection.

Using the Ferno Mondial Funeral Transporter significantly reduces physical strain by transferring the weight of the deceased onto the trolley rather than the caregiver's back. You can wheel the trolley—along with the stretcher and the deceased—directly to the vehicle and effortlessly load the entire unit. The trolley is designed for a loading height of approximately 70 cm.

**Enhanced Maneuverability and Comfort:** A standout feature of the Ferno Mondial Funeral Transporter is the swivel function for all four wheels. By simply activating a lever at the foot end, the front castor wheels are unlocked—enabling full 360° maneuverability. When the front wheels are aligned straight again, they automatically lock back into place. This allows you to easily and precisely position the trolley close to beds or narrow spaces.

**3 in 1 solution:** Thanks to its high load capacity, the Ferno Mondial Funeral Transporter is not limited to use with stretchers. It is also suitable for coffin transport. With the included multifunctional adapter, both a stretcher and a coffin can be securely positioned and transported on the same stretcher, making a truly versatile solution for funeral homes and mortuary services.

- Color-coded controls with intuitive operation and automatic leg release
- Large wheels for smooth rolling on all surfaces
- Front wheel locking system for optimal steering control in confined spaces
- Six adjustable height levels for ergonomic handling and easier transfer of the deceased
- Durable metal frame for long-lasting, high-capacity use
- Ergonomically shaped controls at both the head and foot ends for comfortable handling
- Two loading rollers at the head end for easier vehicle access

Ferno Mondial Funeral Transporter is complete with adapter F0M161S and the M3369 stop plate.

Recommended option: F011U4 - Mortuary Stretcher M11U foldable with grey cover

<b>Art. nr</b>	F0M6F-M3369
<b>Height</b>	35-103 cm
<b>Length</b>	203 cm
<b>Load Capacity</b>	270 kg
<b>Weight</b>	36 kg
<b>Width</b>	59 cm



## Mortuary Cot M11U Foldable with Grey Cover

FERNO transfer stretchers are made from high-quality aluminum, making them lightweight yet highly durable. The stretcher comes with two foldable support legs, two foldable wheels, a foot pocket, and a footrest as standard. The gliding rollers at the foot end allow for easy maneuverability, even when the stretcher is raised, such as on stair landings in narrow staircases.

The transfer stretcher can be folded in the middle, allowing for space-saving storage. The base covering is washable and resistant to disinfectants. The body cover varies depending on the model. All stretcher models are loadable up to 160 kg.

Model 11 U features a black base covering. The transport cover is attached to the stretcher frame using Velcro and six holding loops, making it easy to remove or replace. Cut-outs for the three standard straps are integrated. The transport cover is plastic-coated on the inside and available in four different colors on the outside.

<b>Art. nr</b>	F011U4
<b>Height</b>	21 cm (8 cm folded)
<b>Length</b>	190 cm (103 cm folded)
<b>Load Capacity</b>	160 kg
<b>Weight</b>	11 kg
<b>Width</b>	54 cm



## Mortuary Cot M11BS Wide Body


Designed specifically for the transport of very heavy individuals, the new, wider, and longer Transfer Stretcher Model 11 BS provides a reliable solution for handling large patients.

Like the proven FERNO transfer stretchers, this model is made from a robust aluminum frame. A metal plate increases stability for added support.

With a width of 60 cm, it can easily accommodate overweight patients. The standard features include two short support legs, gliding rollers at the foot end, and a foldable footrest for added convenience.

While the stretcher can be stored flat for space-saving purposes, it is not foldable for stability reasons. The package includes three securing straps for safe transport.

<b>Art. nr</b>	F011BS
<b>Height</b>	10 cm (with footrest folded: 15 cm, with footrest raised: 24 cm)
<b>Length</b>	199 cm
<b>Load Capacity</b>	250 kg
<b>Weight</b>	10 kg (12 kg with cover and straps)
<b>Width</b>	60 cm



# MORTUARY STRETCHERS

## Mortuary Cot M24MM Heavy Duty w. Wheel Locks

First call stretcher for easy handling of the heavy lifts. A new and innovative design ensures ergonomic and safe handling of even the heaviest loads. Simple and safe loading into a vehicle.

- Independent height adjustment in head-and foot end.
- Multi level height adjustment.
- Stands up for the heaviest loads.
- Big wheels for ease of maneuverability.
- Supplied with 2 belst and a mattress.
- A very compact and solid cot.
- Black powder paint of the undercarriage.
- Loading height 61-66 cm.

<b>Art. nr</b>	F024MM
<b>Height</b>	99 cm
<b>Length</b>	196 cm
<b>Load Capacity</b>	450 kg
<b>Weight</b>	33 kg
<b>Width</b>	58 cm

## Mortuary Cot M24 with Sidearms & Wheel Locks

First call stretcher. Simple to use and ergonomic to work with. Many high positions. Head- and foot end to be operated seperately avoiding heavy lifts. Delivered with a thin mattress and belts.

- Can be supplied in an ekstra heigh version for up to 78 cm loading height.
- Sidearms can be supplied as an option when ordering a new cot.
- A special transfer mattress with lifting handles can be used on this cot (optional).

<b>Art. nr</b>	F024
<b>Height</b>	82 cm (24 cm folded)
<b>Length</b>	196 cm
<b>Load Capacity</b>	182 kg
<b>Weight</b>	24,5 kg
<b>Width</b>	52 cm

## 11T Funeral Track Stretcher

Foldable stretcher for easy handling on stairs. A foldable belt system is fixed to the stretcher. The belt system is pulled out and seriously eases the sliding down stairs.

The system folds into a compact unit.

Supplied with 2 two part belts.



<b>Art. nr</b>	F011T
<b>Height</b>	19 cm (15 cm folded)
<b>Length</b>	187 cm
<b>Load Capacity</b>	159 kg
<b>Weight</b>	10,6 kg
<b>Width</b>	47 cm

## Funeral Stretcher Foldable Black

The Model 11 Funeral Stretcher is designed to streamline transportation and storage while maintaining the dignity and professionalism required in funeral services. Its foldable design allows for easy handling in a variety of environments, from streets to vehicles, ensuring a seamless workflow.

Lightweight yet durable, this stretcher provides the reliability you need during critical moments.

- **Effortless Maneuverability:** Equipped with folding posts and wheels for smooth handling and easy navigation in tight spaces.
- **Compact and Convenient Storage:** Slimline hinges enable the stretcher to fold neatly in half, optimizing storage space in vehicles or facilities.
- **Secure and Reliable:** Includes two 5-foot quick-release restraints to ensure secure transport.
- **Durable Construction:** Made with a lightweight aluminum frame and a sturdy, vinyl-coated nylon cover, built to withstand rigorous use while being easy to clean.

<b>Art. nr</b>	F011BLACK
<b>Height</b>	47 cm
<b>Length</b>	187 cm
<b>Weight</b>	7 kg
<b>Width</b>	47 cm



# MORTUARY TABLES

## Mortuary Hydraulic Operating/Transfer Table M101H SS

Allows quick height and angle adjustments using the three foot pedals. Either end of the table may be raised or lowered to any position between 28½" and 40½" above the floor for easy preparation from a standing or sitting position. The table top is convex to drain fluids away from the body and an extra-deep trough runs the perimeter of the table, reducing the need for splash guards.

<b>Art. nr</b>	F0661010
<b>Height</b>	102.9 cm
<b>Length</b>	203.2 cm
<b>Load Capacity</b>	272 kg
<b>Weight</b>	72 kg
<b>Width</b>	76.2 cm

## Mortuary Foldable Operating/Transfer Table M102 SS

A true combination table, ideal for embalming, dressing and transferring. An exclusive, one-hand release easily adjusts either end of the table from 32½" to 38½" with the ends being independently adjustable. The table is constructed of a stainless steel top with a convex surface, tubular aluminum frame, and 6" wheels with ball bearings.

<b>Art. nr</b>	F0661020
<b>Height</b>	98 cm
<b>Length</b>	205 cm
<b>Load Capacity</b>	454 kg
<b>Weight</b>	40 kg
<b>Width</b>	77 cm

## Mortuary Operating/Transfer Table M103 SS

A multipurpose table, ideal for embalming, dressing and transferring. For added convenience, each end can be adjusted independently for individual working height and drainage. The stainless steel top is durable and rigid, featuring a convex surface and deep side troughs for efficient drainage. Six-inch wheels with locks keep the table stationary while in use.

<b>Art. nr</b>	F0661030
<b>Height</b>	95 cm
<b>Length</b>	205 cm
<b>Load Capacity</b>	454 kg
<b>Weight</b>	36 kg
<b>Width</b>	77 cm



## Folding Dressing Table M34 Gray (Regular/Wide)

The Model 34 Dressing Table unfolds quickly for use and stores compactly when not in use. The two-piece, laminated top resists scratches and stains and can be cleaned in seconds. Four swivel wheels are equipped with locks.

Folds for storage, easy-to-clean table top and six-inch swivel wheels.

<b>Art. nr</b>	F034 (Wide: F0663920)
<b>Height</b>	155 cm
<b>Length</b>	196 cm
<b>Load Capacity</b>	227 kg
<b>Weight</b>	31 kg
<b>Width</b>	62 cm (Wide: 74 cm)



## Folding Dressing Table M36 Gray (Regular/Wide)

An exclusive, one-hand release on the Model 36 allows you to quickly adjust the table height from 81 cm to 97 cm. The Model 36 rolls on four swivel wheels that lock directionally for stability.

Adjustable Height, easy-to-clean table top and six-inch swivel wheels.

<b>Art. nr</b>	F0663602 (Wide: F0663924)
<b>Height</b>	97 cm
<b>Length</b>	196 cm
<b>Load Capacity</b>	227 kg
<b>Weight</b>	35 kg
<b>Width</b>	62 cm (Wide: 74 cm)



## Mortuary Cot, Mini 141cm

Compact and lightweight, designed for pet removal or small body retrieval. Its One-Man®, roll-in style allows one person to perform retrieval and has extendable loading wheels for easy loading and unloading. Special release levers make it possible to raise and lower the cot without assuming the full weight of the cot and body.

<b>141 cm</b>	F0403026
<b>Height</b>	75 cm
<b>Length</b>	141 cm
<b>Load Capacity</b>	136 kg
<b>Weight</b>	28.5 kg
<b>Width</b>	56 cm





Mortuary Cemetery Church Truck M87CCT, Silver w. Handles

The Ferno 87 Cemetery Church Truck is designed for graveside use, with large pneumatic wheels, durable frame, and rubber bumpers for casket security. This church truck comes with four side handles and two wheel locks on all configurations.

The Ferno 87 Cemetery Church Truck can easily overcome grass and graveside terrain smoothly and reliably, thanks to the aluminum frame and ridged wheels.

Two of the 8-inch wheels are equipped with wheel locks for stability and the undercarriage itself can be adjusted to four different positions to accommodate caskets of different size. Without a casket mounted, this church truck weighs just under 60 lbs, and can be quickly transported around a complex.

**Length:** Folded: 17 cm, 1st Position: 50 cm, 2nd Position: 75 cm, 3rd Position: 128 cm, Fully Extended: 156 cm  
**Height:** Folded: 83 cm, Fully Extended: 64 cm

Art. nr	F087HCCT
Load Capacity	454 kg
Weight	27 kg
Width	46 cm



Mortuary Cemetery Church Truck M87CCT, Black w. Handles

The Ferno 87 Cemetery Church Truck is designed for graveside use, with large pneumatic wheels, durable frame, and rubber bumpers for casket security. This church truck comes with four side handles and two wheel locks on all configurations.

The Ferno 87 Cemetery Church Truck can easily overcome grass and graveside terrain smoothly and reliably, thanks to the aluminum frame and ridged wheels.

Two of the 8-inch wheels are equipped with wheel locks for stability and the undercarriage itself can be adjusted to four different positions to accommodate caskets of different size. Without a casket mounted, this church truck weighs just under 60 lbs, and can be quickly transported around a complex.

**Length:** Folded: 17 cm, 1st Position: 50 cm, 2nd Position: 75 cm, 3rd Position: 128 cm, Fully Extended: 156 cm  
**Height:** Folded: 83 cm, Fully Extended: 64 cm

Art. nr	F0618887
Load Capacity	454 kg
Weight	27 kg
Width	46 cm

## Mortuary Church Truck M87H, Silver w. Handles

Made of strong, rectangular, aluminum tubing, the Model 87 Church Truck features a protective anodized finish that keeps the truck looking bright and shiny-new as well as making it easy to clean.

The stable x-frame undercarriage and rubber rests help to secure the casket. Six-inch rubber tire roll easily over carpet, sills and other obstacles with quiet dignity. When moved without a casket a special built-in handle allows the truck to be pushed, pulled or carried.

The Model 87 Church Truck locks into four positions, allowing a wide range of casket sizes to be displayed with one truck. The additional length also increases the stability and makes an attractive appearance.

The Model 87-H Church Truck is identical to the Model 87 but with the added convenience of fold-out carrying handles.

<b>Art. nr</b>	F087H
<b>Height</b>	53 cm (73 cm folded)
<b>Length</b>	157 cm (16 cm folded)
<b>Load Capacity</b>	454 kg
<b>Weight</b>	16 kg
<b>Width</b>	55 cm



## Mortuary Church Truck M87H, Black w. Handles


Made of strong, rectangular, aluminum tubing, the Model 87 Church Truck features a protective anodized finish that keeps the truck looking bright and shiny-new as well as making it easy to clean.

The stable x-frame undercarriage and rubber rests help to secure the casket. Six-inch rubber tire roll easily over carpet, sills and other obstacles with quiet dignity. When moved without a casket a special built-in handle allows the truck to be pushed, pulled or carried.

The Model 87 Church Truck locks into four positions, allowing a wide range of casket sizes to be displayed with one truck. The additional length also increases the stability and makes an attractive appearance.

The Model 87-H Church Truck is identical to the Model 87 but with the added convenience of fold-out carrying handles.

<b>Art. nr</b>	F0618869
<b>Height</b>	53 cm (73 cm folded)
<b>Length</b>	157 cm (16 cm folded)
<b>Load Capacity</b>	454 kg
<b>Weight</b>	16 kg
<b>Width</b>	55 cm



## Church Truck Compact Silver M88

Designed for those who prefer a compact truck. The Church Truck is constructed of rectangular aluminum tubing with a protective anodized finish to keep the truck looking bright and shiny new as well as making it easy to clean. This truck features wheels and three locking frame positions.

The Model 88-H Church Truck is identical to the Model 88 but with the added convenience of fold-out carrying handles. Use the handles to navigate sills, curbs and other obstacles or for assistance when climbing stairs.

<b>Art. nr</b>	F088
<b>Load Capacity</b>	454 kg
<b>Weight</b>	15 kg
<b>Width</b>	56 cm
<b>Height</b>	55 cm (74 cm folded)
<b>Length</b>	121 cm (19 cm folded)

## Church Truck Compact Silver w/Handle M88

Designed for those who prefer a compact truck. The Church Truck is constructed of rectangular aluminum tubing with a protective anodized finish to keep the truck looking bright and shiny new as well as making it easy to clean. This truck features wheels and three locking frame positions.

The Model 88-H Church Truck is identical to the Model 88 but with the added convenience of fold-out carrying handles. Use the handles to navigate sills, curbs and other obstacles or for assistance when climbing stairs.

<b>Art. nr</b>	F088H
<b>Load Capacity</b>	454 kg
<b>Weight</b>	15 kg
<b>Width</b>	56 cm
<b>Height</b>	55 cm (74 cm folded)
<b>Length</b>	121 cm (19 cm folded)



## Mortuary Showroom Truck M110 Aluminium

The Model 110 is a dual-duty truck that works as a showroom and utility truck. The polished aluminum frame is exceptionally strong and will look like new for years with minimal maintenance.

Rubber tips hold caskets or rough boxes firmly in place while five-inch ball bearing wheels roll over obstacles with ease.



<b>Art. nr</b>	F110
<b>Height</b>	56 cm
<b>Length</b>	112 cm
<b>Load Capacity</b>	272 kg
<b>Weight</b>	9 kg
<b>Width</b>	51 cm

## Mortuary Showroom Truck XL M112 Aluminium

The extra-long and wide wheelbase on the Model 112 Showroom Truck provides exceptional stability. The Model 112 performs equally well as a utility truck.

The strong aluminum frame will keep its high polish for years with a minimum of maintenance. Rubber tips hold caskets or rough boxes firmly in place and five-inch, ball-bearing wheels roll over obstacles with ease.



<b>Art. nr</b>	F112
<b>Height</b>	54 cm
<b>Length</b>	121 cm
<b>Load Capacity</b>	272 kg
<b>Weight</b>	9 kg
<b>Width</b>	56 cm



## EVACuation Mattress Funeral

The Bodypod® EVACuation Transfer Mattress is a robust, ergonomic, and hygienic solution for the dignified handling of deceased individuals. Designed with both efficiency and user comfort in mind, it is an essential tool for funeral professionals, first responders, and emergency teams who operate in demanding environments. This innovative mattress combines lightweight maneuverability, durability, and ethical considerations.

<b>Art. nr</b>	5006056
<b>Height</b>	7 cm
<b>Length</b>	210 cm
<b>Load Capacity</b>	200 kg
<b>Weight</b>	6 kg
<b>Width</b>	60 cm



## Ferno Scoop Stretcher 65EXL

The Scoop Stretcher is designed for the dignified and efficient lifting and transfer of deceased individuals. The 65EXL model stands out for its high level of safety and stability, lightweight construction, ease of cleaning, and user-friendly design.

Simply scoop up the deceased and transfer them to a transport stretcher or coffin. Three two-piece restraint straps are included in the standard delivery.

<b>Art. nr</b>	F065EXLB
<b>Height</b>	8 cm
<b>Length</b>	165 - 201 cm
<b>Load Capacity</b>	227 kg
<b>Weight</b>	8 kg
<b>Width</b>	43 cm



## Germa Easy Handle 4-kit

Germa Easy Handle, the simple solution that makes the big difference. Easy Handle is an ergonomic aid intended for use when lifting off mattresses, mini stretchers, vacuum mattresses, scoop stretchers etc.

Easy Handle enables more ergonomic lifting from knee height. Ferno Easy handle kit has 4 slings with ergonomic handles. Everything packed in a practical bag.

<b>Art. nr</b>	GE15000000034
<b>Load Capacity</b>	250 kg
<b>Weight</b>	0,45 kg
<b>Size</b>	Bag 20x30 cm

## Funeral Transfer Mattress Black

Transfer Mattress Black is designed specifically for the funeral market as an evolution of the classic acute mattress used in the ambulance services. Made of high quality material, is compatible with most funeral stretchers on the market, has 6 pieces of ergonomic handles on each side of the mattress, makes the mattress easy to handle and to lift and carry a body.

<b>Art. nr</b>	GE271008800
<b>Height</b>	3,5 cm
<b>Length</b>	195 cm
<b>Load Capacity</b>	250 kg
<b>Weight</b>	3,8 kg
<b>Width</b>	50 cm



## Step by Step Black

Step by Step is quick and easy to apply to the mattress and includes a restraint system for the patient to be secured on the mattress.

This solution delivers a secure evacuation of laying patients with ergonomics and less hassle for the rescue teams.

<b>Art. nr</b>	GE15000000117
<b>Height</b>	15 cm
<b>Length</b>	30 cm
<b>Load Capacity</b>	250 kg
<b>Weight</b>	1,3 kg
<b>Width</b>	20 cm



## Germa Minitranz, Mini stretcher

To make usage as easy as possible, the back of the Minitranz mattress is made of friction-reducing material to slide easily behind or under the patient. Glass fibre material on the back is very durable and can be used sliding on floor or ground. The great benefit is that you are able to bend the flexible stretcher/mattress around the patient and in step two lift and move the patient carrying in 10 possible strong handles.

<b>Art. nr</b>	GE70050
<b>Height</b>	1,5 cm
<b>Length</b>	136 cm
<b>Load Capacity</b>	250 kg
<b>Weight</b>	1,8 kg
<b>Width</b>	76 cm (46 cm folded)





## Lifting Handle for Scoops

The carrying straps minimize the need to bend when lifting scoop stretchers, backboards, and basket stretchers.

The carrying straps are simply looped around the handle recesses. The integrated ergonomic handle makes lifting and carrying easier.

<b>Art. nr</b>	F0050135
----------------	----------



## Ferno Transfer Board 114 cm

Use the Ferno® Transfer Board to transfer a patient between adjoining surfaces, such as a bed and a transport cot, or a transport cot and a table. High density polyethylene construction makes the transfer board light and durable. Six handles ensure that you always find a good grip. X-ray transparent. Well suited for the larger patients.

<b>Art. nr</b>	F673
----------------	------

<b>Height</b>	0,6 cm
---------------	--------

<b>Length</b>	114 cm
---------------	--------

<b>Weight</b>	2 kg
---------------	------

<b>Width</b>	46 cm
--------------	-------



## Glide Board 150 kg

The Glide Board is designed for the safe and efficient repositioning of individuals with a maximum weight of 150 kg. The tubular outer material features an internal high-sliding specialty fabric that enables smooth movement. A 3 mm thick plastic plate inserted into the tube supports the rolling process. This design allows the repositioning of a person without the need to lift them.

<b>Art. nr</b>	FORBF
----------------	-------

<b>Height</b>	0.5 cm
---------------	--------

<b>Length</b>	75 cm
---------------	-------

<b>Load Capacity</b>	150 kg
----------------------	--------

<b>Weight</b>	1.16 kg
---------------	---------

<b>Width</b>	45 cm
--------------	-------



## Mortuary Cot Mattress 178x55x3cm for M24MM Black

These mattresses are constructed with a one-inch polyurethane, foam pad cover with vinyl-coated nylon.

The cover is wear and abrasion-resistant, bacteriostatic, fire-resistant, and easy to clean.

Art. nr F0314042



## Funeral Black Mattress for M11BS

Custom-Fitted Black Mattress for Model 11BS Stretcher. Designed specifically for the Model 11BS mortuary stretcher, this black mattress offers a perfect blend of functionality and professionalism. Crafted from a highly durable, easy-to-clean material, it ensures maximum hygiene even in the most demanding environments.

The form-fitted design guarantees a secure, tailored fit that not only enhances the stretcher's appearance but also provides optimal protection and support during transport. Built to withstand daily use, this premium mattress is an essential addition for any funeral service provider looking to maintain high standards of care, dignity, and presentation.

Art. nr F0359BS



## Germa Easy Lift Shoulder

Easy Lift Shoulder is a newly developed tool that relieves the hands and arms when lifting or carrying a stretcher. Easy Lift is easy to use and has unique features that make it an efficient ergonomic tool for the user.

Easy Lift has friction lock and is providing a firm grip with very limited power required from the user.

Art. nr GE15000000016

Weight 0,4 kg



## Mortuary Church Truck Drape Frame M87F

The Model 87-F Drape Frame Attachment is a metal lattice frame that attaches to the church truck at six points to provide strong, yet lightweight drape mounting.

Art. nr F0648700

## Mortuary Church Truck Drape For M87 Crushed Velvet Red

The Model 87-D Velvet Drapes fit the Model 87 Church Trucks.

Art. nr F0312890

## Mortuary Church Truck Drape For M87 Crushed Velvet Black

The Model 87-D Velvet Drapes fit the Model 87 Church Trucks.

Art. nr F0314032

## Mortuary Church Truck Drape For M87 Crushed Velvet Blue

The Model 87-D Velvet Drapes fit the Model 87 Church Trucks.

Art. nr F0312900



## Church Truck Cover Crypton® Aria fabric M402 W.Drapes Red

The Ferno church truck cover allows you to protect your church truck from dust and scratches while in storage. The Model 402 is designed for use on church trucks with drapes.

Model 402 dimensions: 66 x 70 x 35 cm

Art. nr F0314150



## Church Truck Cover Crypton® Aria fabric M402 W.Drapes Black

The Ferno church truck cover allows you to protect your church truck from dust and scratches while in storage.

The Model 402 is designed for use on church trucks with drapes.

Model 402 dimensions: 66 x 70 x 35 cm

Art. nr F0314153



# DRAPES AND COVERS



## Cot Pouch Suede Zipper 191x58cm M335 Red/Black

The First Call Pouch comes complete with a bottom and full-length brass zipper.

The pouch is secured by snap fastener tabs and features full expansion restraint strap slots. Fits One-Man® cots and stretchers approximately 191 cm x 58 cm.

**Art. nr** F0314075 (red), F0314076 (black)



## Cot Cover Slip over Suede 193x51cm M320 Red/Black

The slip-over styling of this cover fits One-Man® Cots. Secured by elastic at both ends, the cover is reversible end-to-end.

The cover fits cots approximately 76" x 20" (193 cm x 51 cm), which includes the Model 24 series and Model 27-1 series cots.

It is available in the Crypton® Aria or Suede. Includes a 12" x 14" pocket.

**Art. nr** F0314081 (red), F0314082 (black)



## Cot Cover Slip over Suede 193x58cm M321 Red/Black

This cover features tailored, slip-over styling and is reversible end-to-end.

Fits cots 76" x 23" (193 cm x 58 cm), which includes the Model 24 series and Model 23 series cots.

It is available in the Crypton® Aria or Suede. Includes a 12" x 14" pocket.

**Art. nr** F0314087 (red), F0314088 (black)



## Cot Cover Suede Zipper 191x58cm M330 Red/Black

This cover fits One-Man® cots, most stretchers and standard cots and has a full-length brass zipper.

The cover features snap-fastener tabs that hold it in place beneath the mattress and it has full-expansion restraint strap slots.

Fits cots approximately 75" x 23" (191 cm x 58 cm).

It is available in the Crypton® Aria or Suede. Includes a 12" x 14" pocket.

**Art. nr** F0314093 (red), F0314094 (black)



## Cot Cover Slip over Crypton® Aria 193x51cm M320 Red/Black

The slip-over styling of this cover fits One-Man® Cots. Secured by elastic at both ends, the cover is reversible end-to-end

The cover fits cots approximately 76" x 20" (193 cm x 51 cm), which includes the Model 24 series and Model 27-1 series cots.

It is available in the Crypton® Aria or Suede. Includes a 12" x 14" pocket.

**Art. nr** F0314126 (red), F0314129 (black)



## Cot Cover Slip over Crypton® Aria 193x58cm M321 Red/Black

This cover features tailored, slip-over styling and is reversible end-to-end.

Fits cots 76" x 23" (193 cm x 58 cm), which includes the Model 24 series and Model 23 series cots.

It is available in the Crypton® Aria or Suede. Includes a 12" x 14" pocket.

**Art. nr** F0314130 (red), F0314133 (black)



# DRAPES AND COVERS



## Cot Cover Crypton® Aria Zipper 191x58cm M330 Red/Black

.....

This cover fits One-Man® cots, most stretchers and standard cots and has a full-length brass zipper.

The cover features snap-fastener tabs that hold it in place beneath the mattress and it has full-expansion restraint strap slots.

Fits cots approximately 75" x 23" (191 cm x 58 cm).

It is available in the Crypton® Aria or Suede. Includes a 12" x 14" pocket.

**Art. nr** F0314138 (red), F0314141 (black)



## Cot Pouch Crypton® Aria Zipper 191x58cm M335 Red/Black

.....

The First Call Pouch comes complete with a bottom and full-length brass zipper.

The pouch is secured by snap fastener tabs and features full expansion restraint strap slots.

Fits One-Man® cots and stretchers approximately 191 cm x 58 cm.

**Art. nr** F0314142 (red), F0314145 (black)



## Mortuary Cot Cover Crypton® Aria fabric for Mini, Dark red

.....

Cover designed specifically for the Mini Cot Roll-In Style Cot.

**Art. nr** F0314158

## Cover for 11U Black 190x51 cm

.....

This cover is available in the current length of 190 cm as standard. It is smooth with sewn-on Velcro strips and comes in black. The cover does not include protective covers and is designed for the attachment of transport covers (which are fixed using the sewn-on Velcro).

The cover is also designed to be attached to the stretcher in the same way as the Model 11.

The image shows the stretcher with the black cover and one of the transport covers described on the following page. Please note that the transport cover is not included with the replacement cover.

Art. nr F0EB115



## Cover for 11U Blue 190x51 cm

.....

This cover is available in the current length of 190 cm as standard. It is smooth with sewn-on Velcro strips and comes in black. The cover does not include protective covers and is designed for the attachment of transport covers (which are fixed using the sewn-on Velcro).

The cover is also designed to be attached to the stretcher in the same way as the Model 11.

The image shows the stretcher with the black cover and one of the transport covers described on the following page. Please note that the transport cover is not included with the replacement cover.

Art. nr F0335U1



## Cover for 11BS Blue 190x61 cm

.....

The transport cover 335 UB has the properties of the transport cover 335 U, but is dimensionally adapted to the transfer stretcher 11 BS.

Art. nr F0335UB1



## Body bag

The VHN 20 body wraps are made of opaque, white, 0.15 mm thick plastic film and feature a nylon zipper.

Size approx.: 220 x 90 cm

<b>Art. nr</b>	F0VHN20
----------------	---------

## Body Bag Disposable PEVA Burnable M347-1 White

This disposable pouch is made of electronically-sealed vinyl and resists abrasions, acid, formaldehyde and blood stains. It will not leak, deteriorate or mildew, and it is not affected by temperature extremes.

Available in white, the Model 347-1 features a heavy-duty, U-shaped metal zipper with double slides.

<b>Art. nr</b>	F3471
----------------	-------

<b>Length</b>	231 cm
---------------	--------

<b>Width</b>	96 cm
--------------	-------

## Body bag XL Absorbent Maroon

Absorbent Extra Large / funeral bag / body bag with an absorbent inside. Made of Polyethylene and equipped with strong handles with a load capacity of 180kg.

Environmentally adapted, which means that the bag can be buried in its entirety and also used for incineration because the plastic is converted into water vapor. Shelf life for several years as long as it is not exposed to direct sunlight.

<b>Art. nr</b>	5001009
----------------	---------

<b>Length</b>	255 cm
---------------	--------

<b>Load Capacity</b>	180 kg
----------------------	--------

<b>Weight</b>	1,7 kg
---------------	--------

<b>Width</b>	114 cm
--------------	--------



## Body bag Absorbent Standard Maroon

Absorbent standard body bag / funeral bag / body bag with an absorbent inside. Made of Polyethylene and equipped with strong handles with a load capacity of 180kg.

Environmentally adapted, which means that the bag can be buried in its entirety and also used for incineration because the plastic is converted into water vapor. Shelf life for several years as long as it is not exposed to direct sunlight.

<b>Art. nr</b>	5001001
<b>Length</b>	215 cm
<b>Load Capacity</b>	180 kg
<b>Weight</b>	1,5 kg
<b>Width</b>	95 cm



## Body bag Standard Maroon

Absorbent Extra Large / funeral bag / body bag with an absorbent inside. Made of Polyethylene and equipped with strong handles with a load capacity of 180kg.

Environmentally adapted, which means that the bag can be buried in its entirety and also used for incineration because the plastic is converted into water vapor. Shelf life for several years as long as it is not exposed to direct sunlight.

<b>Art. nr</b>	5001006
<b>Length</b>	215 cm
<b>Load Capacity</b>	180 kg
<b>Weight</b>	1 kg
<b>Width</b>	95 cm



## Body bag XL Maroon

Standard body bag / funeral bag / body bag without absorbent inside. Made of Polyethylene and equipped with strong handles with a load capacity of 180kg.

Environmentally adapted, which means that the bag can be buried in its entirety and also used for incineration because the plastic is converted into water vapor. Shelf life for several years as long as it is not exposed to direct sunlight.

<b>Art. nr</b>	5001008
<b>Length</b>	255 cm
<b>Load Capacity</b>	180 kg
<b>Weight</b>	1,4 kg
<b>Width</b>	114 cm



## OTHER MORTUARY PRODUCTS

### Mondial Funeral Lock

Floor-Mounted Lock for Ferno Mondial Funeral Transporter.

This durable floor-mounted lock is specifically designed to secure stretchers during transport in Hearses. Built with strength and reliability in mind, the lock ensures that the stretcher remains firmly in place.

Its robust construction and wear-resistant materials offer long-lasting performance, while the simple locking mechanism makes it quick and easy to operate. Discreet in appearance and user-friendly in design, this lock meets the high standards of the funeral industry for safety, hygiene, and professionalism.

**Art. nr** F0FF50SK

### Mortuary Fastener M5 for funeral stretchers

Permanently mounted stretcher fastener to secure cots and stretchers. Compatible with the stretchers model 24, 24MAXX and 24MinMaxx.

Thumbscrews permit easy installation and removal.

The rear fastener hook stores flat and rotates up to function as a ratchet mechanism when needed.

Not for use in crash stable applications.

**Art. nr** F005

### Side Lift-Handles (pair, clamp-on) For M24MAX/ M24MM

The Ferno Swing Out Lifting Handles can be attached to the mainframe of the cot, to provide more lifting points for emergency responders to grasp and lift the cot.

This is especially useful for taller or larger patients.

**Art. nr** F0821998

### Wheel Cups M533 (each)

Designed to stop stretchers moving around inside a vehicle.

Polished aluminum.

Available in 4" (5 cm) , 5" (7,5cm) and 7" 12,5 cm)

**Art. nr** F533

**Material** Aluminium



### Recessed Post Cups M534 (pair)

A secure way to hold stretcher posts in place in the vehicle floor.

Permanently mounted recessed in the vehicle.

**Art. nr** F534

**Package size** 2 pcs



### Raised Post Cups M535 (pair)

A secure way to hold stretcher posts in place in the vehicle floor.

Permanently mounted recessed in the vehicle.

**Art. nr** F534

**Package size** 2 pcs





**Ferno Norden**, founded in 1979, is a leading provider of specialized equipment for prehospital care, industrial safety, rescue, and funeral services across the Nordic region. With decades of experience and a commitment to safety, quality, and innovation, our products are trusted by emergency services, defense forces, industrial operators, and healthcare providers. Our extensive offerings also include advisory services and customized installations for specialized vehicles.



With dedicated teams and offices across several European countries, Ferno Norden is committed to meeting the safety, rescue, and emergency needs of the Nordic and Baltic regions. Ferno Norden owns and operates its sales offices in **Norway**, **Sweden**, **Denmark**, and **Finland**, while trusted local sales teams represent us in the other countries such as **Iceland**, **Lithuania**, **Latvia** and **Estonia**, ensuring consistent expertise and support tailored to each region. For customers in the **Faroe Islands** and **Greenland**, sales and support are handled through our Danish office.



Whether you need product information, assistance with an order, or guidance on selecting the right solutions for your organization, our teams are ready to help. Contact your local office to find the best tools and systems to support your mission.

[www.fernonorden.com](http://www.fernonorden.com)

## *Our sustainability work*

Our sustainability work follows and is based on the UN's Agenda 2030, which consists of 17 global goals for sustainable development that aim to eradicate poverty, stop climate change and create peaceful and safe societies.

