



CHE UNITED KINGDOM VEHICLE APPROVAL AUTHORITY

**BS EN 1789:2007 +A2 1789:2014  
Certificate of Compliance**

This component has satisfied all requirements of paragraph(s) 4.5.9 and 5.4 of BS EN 1789:2007 including Amendment +A2 1789:2014 in order to achieve CEN compliance for Dynamic testing. Appropriate tests have been witnessed by the Vehicle Certification Agency and the specification of the component tested has been documented.

See VCA test report: EST402492

**Reason for Extension:** Not applicable

- 1. Make (trade name of manufacturer):** Ferno
- 2. Type and general description:** Pedimate Plus; Pedimate Plus Child Harness
- 3. Means of Identification:** Labelling and Serial Number
- 4. Location of that Marking:** Rear side
- 5. Category of Vehicle:** Not applicable
- 6. Name and Address of Manufacturer:**  
Ferno UK Limited  
c/o Ferno-Washington, Inc.  
70 Weil Way  
Wilmington  
OH 45177  
United States
- 7. If approved as a Separate Technical Unit, Location of that marking:** Rear side
- 8. Remarks:**
  - a) An information file containing test reports, drawings and statements is held on file at VCA.
  - b) Tested on PowerX Stretcher top with loops on harness fixations (triangular bracket locating into stud); back rest at 45 degrees. Nominal max load mass
  - c) Results may be read across to the standard Pedimate and Neomate devices

**Date:** 21 MARCH 2022

**C MCCABE**  
Chief Technical and Statutory Operations Officer

EST402492

