

PRODUKTINFORMATION
STOP THE BLEED KIT SCOOP AND RUN COMPLETE BAG 10 SET
Productno.: 5020208



STOP THE BLEED KIT SCOOP AND RUN COMPLETE BAG 10 SET

A complete kit for mass casualty's packed in a well-organized bag.

Designed and packed for grab and go situations

The main space is filled with 10 sets of vacuum-packed Stop the Bleed kits including:

1 CAT Tourniquet, 1 blanket, 1 compressed gauze and 1 trauma bandage.

The bag also includes 3 evacuation stretchers and 2 splints

Product information:

The bag has several pockets and 6 Velcro attached inner bags. The design is compact, and you can stack them on each other without any problems. Handles on top and a shoulder harness. Rubber reinforced plates to protect the bag from the ground. Bright red color and reflexes.

Certifications:

CE

Attributes

Length: 45cm

Width: 35cm

Height: 23cm

Weight: 6 kg packed

Color: Red