

# NORSE RESCUE® Eye Shield

## Red - One Size

---



---

Item No:	CL-9
Manufacturer:	Gulden Ophthalmics
Certificates:	CE
Product Size:	7,3 x 6 x 0,3 cm
Color:	Red
Country of origin:	US
NSN no:	6515-01-647-5456
Tariff:	90019070
Weight:	7,5 g

The NORSE RESCUE® Eye Shield is a Polycarbonate Eye Shield that meets current military pre-hospital guidelines for the treatment of traumatic eye injury in the field. This convex rigid eye shield allows the user to protect the injured eye and blocks external pressure from being applied directly to the injury, preventing further damage.

Made of a high strength, translucent red polymer, it is shatterproof and will not flatten or disfigure when pressure is applied. Holes in the eye shield provide ventilation, and is designed for use over either eye.