

NORSE RESCUE® Splint

Flat Folded



Item No:	IM-00047
Manufacturer:	EMSRUN
Certificates:	CE
Package Size:	16 x 11 x 3 cm
Product Size:	36" / 91,4cm
Color:	Black
Country of origin:	Denmark/China
NSN no:	N/A
Tariff:	30059099
Weight:	150g

This NORSE RESCUE® Splint is a aluminum splint suitable for splinting all extremities. It is made of lightweight and malleable foam padded aluminum, it is soft and adaptable. The splint becomes rigid and supportive for fractured or injured limbs by shaping it into a simple C, W or T-shaped curve.

Can be cut with ordinary household scissors to fit any size.

Details:

- Can be rolled or folded for easy storage
- Waterproof
- Lightweight and compact
- Reusable