

Ferno Power X2 Ambulance stretcher

POWERX2™



PRODUCT INFO

Product Nr:	F14000
Height:	
Length:	165 cm / 200 cm
Width:	
Weight:	
Max Load:	318 kg / 498 kg
Compliance:	EN 1789 / EN 1865

The **Ferno Power X2 Ambulance Stretcher** sets a new standard in the industry with its cross-compatibility, enhanced power, improved visibility, and superior maneuverability—outperforming today’s leading powered stretchers. Designed, engineered, and built with safety and ease of use in mind, it excels in even the most challenging circumstances. The **PX2** stretcher provides exceptional stability in unpredictable environments, enabling operators to focus on delivering safe, effective patient care and transportation. With its innovative design, the Ferno **PX2** elevates safety, performance, and maneuverability to new heights. Prioritizing control and reliability, this advanced stretcher empowers healthcare professionals to focus on what truly matters: patient well-being. Resulting from years of meticulous development, the **Ferno Power X2 Ambulance Stretcher** not only meets the demands of the future but exceeds them. It delivers a high-end solution that surpasses today’s standards, ensuring optimal performance and reliability in any situation.



- **CROSSLOCK:** UNIVERSALLY CROSS-COMPATIBLE WITH THE MAJORITY OF MANUAL AND POWER-LOAD SYSTEMS
- **SCENE BRIGHT:** 360° ILLUMINATION, TWO-STAGE LIGHTING
- **EZGLIDE:** 27% LARGER WHEELS, INCREASED GROUND CLEARANCE, ENHANCED STABILITY, GREATER MANEUVERABILITY
- **ENDURACHARGE:** 36-VOLT LiFePO4 BATTERY, LASTS 2X LONGER THAN COMPETITOR'S

STANDARD FEATURES

- Fold Forward Side Arms
- Color Coded Restraints
- Mattress
- Equipment Hook
- Foot-end Wheel Locks
- Foot-end Pull Handle
- Scene Lights
- Status Display with Cycle Counter
- 36V Lithium Battery

AVAILABLE OPTIONS

- 02 Frame Mount
- SX Surface Extenders
- Directional Wheel Locks
- Knee Gatch
- Extending Backrest
- Biosafe Restraints
- IV Pole & Pac-Rac mounts
- Pusch bar
- Storage bag / Storage net