



Ferno F2 Monobloc H56 incl mattress & belts



PRODUCT INFO

Product Nr:	F600341005SN
Length:	193,5 cm
Width:	56,5cm
Height:	28-81cm
Loading:	56 -59 cm
Weight:	43 kg
Max Load:	250 kg
Compliance:	Class 1 MDR / CE
Compliance:	EN1789 & EN1865
Color:	Grey & White

The Ferno F2 Monobloc provides a safe, ergonomic, and reliable patient handling solution. Made of lightweight, durable aluminum, the F2 Monobloc features an independent leg design that supports the patient during loading and unloading, reducing strain on medics. Automatic locking front and rear legs automatically unfold. The color-coded left- and right-hand control levers are user-friendly. Four three-position, telescoping handles provide ergonomic carrying and lifting options for medics. Adjustable six-position backrest provides support for patient comfort, and optional swing-down side rails provide patients a feeling of security.

Option: Ferno F2 Fastener 10G - F600150003

STANDARD FEATURES

- **Excellent maneuverability over any terrain thanks to large 200 mm wheels, control-end swivel lock wheels, and loading-end fixed wheels with wheel locks**
- **Adjustable multi-position backrest supports a variety of comfort positions**
- **Fold-downside rails give patients a feeling of security**
- **Shock frame with contour feature**
- **Low maintenance**
- **Includes mattress, patient restraints, sidearms and a support foot**
- **Telescoping IV-Pole as an option**
- **Color-coded left-hand and right-hand control levers**
- **Loading height in Vehicle 56-59 cm**

