



Ferno Easy Handle 4-kit



Weight	0.45 kg
Color	RED
Load capacity	250 kg
Article number	GE15000000034

Ferno Easy Handle, the simple solution that makes a big difference. Easy Handle is an ergonomic aid intended for use when lifting mattresses, mini stretchers, vacuum mattresses, Scoop stretchers, etc. Easy Handle enables more ergonomic lifting from knee height. Ferno Easy handle kit has 4 slings with ergonomic handles. All packed in a practical bag.

Can withstand a load of 250 kilos.

Washable at 60 degrees for best hygiene.

