



THE BAG - 12 PACK - NR5

The BAG II is a disposable revivator that comes in different sizes.

The BAG II breathing balloon is an excellent choice for simple patient ventilation without the risk of contamination or the need for cleaning. Comes complete with mask and a three meter long oxygen hose. Adjustable handle with hook and loop that fits all types of hands and provides a stable grip. Also includes an oxygen reservoir bag. The children and baby models have a 35 cm H₂O pressure relief valve. Please note that they are sold in packs of 12 pcs / pack.

Art. Nr	L845211
Certifications	MDD/EC/CE
Color	Yellow/Clear
Liquid Capacity	1650 ml
Size	1650 ml
Material	PVC, PC, SI, polypropylene



KEY FEATURES

- Pressure relief valve (35mmH₂O) is available for the models for children and infants
- Does not contain latex
- Dead Space for patient valve: 6.8 ml
- Inhalation resistance: 1.5 cm H₂O at 50 l / min
- Exhalation resistance: 1.8 cm H₂O at 50 l / min
- Temperature limit for storage: -40 ° C to +60 ° C
- Temperature limit for use: -18 ° C to + 50 ° C