



FERNO NATO STRETCHER

NATO Stretcher has been designed to comply with and exceed international standards such as STANAG 2040.

Length (extended handles)	2290 mm
Length (retracted handles)	2000 mm
Width	579 mm
Height	140 mm
Folded size (L/W/H)	230x230x520 mm
Weight	8,5 kg
Load Capacity	250 kg
Supplier Nr	60-0122-001



FERNO NATO STRETCHER

- The most advanced and innovative combat and evacuation system ever designed
- Meets Civilian and Military Fire Safety Codes
- Designed for extreme conditions
- Sturdy Aluminum Frame
- Ergonomically Designed Locking Handles
- 6 IV pole Attachment Points
- 2x Patient-Securing Straps
- Slip-Resistant Pattern Mesh with Fire certification M2/B1, Resistant to low temperatures without cracking (ISO 4675)