



StratLoad™ Attachment Points

The Advantage is an improved redesign from any other plastic style basket stretcher available today. The fully-internal stainless steel frame has a 2500 lb. strength rating while giving a completely smooth sliding surface on the bottom. The outer High Density Polyethylene (HDPE) shell is attached to the top tube, allowing the secondary rail to be completely accessible to provide unlimited options for tying in your patient during the most technical rescue operation.

## Features

- StratLoad™ feature as standard
- Fully-internal stainless steel frame had a 2500 lb. strength rating providing a completely smooth sliding surface
- High Density Polyethylene shell
- 7/8" top tube - The large top tube is very easy to hold onto, especially in long carries
- Closed-Cell Foam Patient Pad -- Will not allow fluids to soak into it
- Recessed sliding surface on the bottom capable of handling rough terrain and protecting the frame
- Class 1 medical device

## Specifications

Advantage Basket Stretcher		
Length	85 in	216 cm
Width	24 in	60 cm
Depth	7.5 in	19 cm
Weight	34.4 lbs	15.6 kg
Rated Load Limit	2500 lbs	11 kN